# Lolina A/S



Address:Sindalsvej 30 8240 Risskov Danmark

Email: Info@lolina.dk

Website:https://lolina.dk

## **Product Specification**

Product name	Lolina® Fibroblast Growth Supplement-2
Cat.No.	NaC102103
Size	5ml
Storage and shipping	Store FGS-2 at $-20^{\circ}$ C; once added to medium, store at $4^{\circ}$ C and protect from light. Do not refreeze after thawing. Dry ice.

## **Product Description**

Lolina® Fibroblast Growth Supplement-2 (FGS-2) is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal fibroblasts. The supplement is formulated (quantitatively and qualitatively) to provide an optimally balanced envi FGS-2 ronment that selectively supports the growth of normal fibroblasts in vitro.

#### **Instructions for Use**

For best results, use FGS-2 in conjunction with 500 ml Fibroblast Medium-2-b and 25 ml fetal bovine serum.

Thaw FGS-2 at 37°C. Gently tilt the tube several times to ensure complete mixing. Spray the medium bottle and supplement tube with 70% ethanol and wipe to remove excess liquid. In a sterile field, remove the caps without touching the interior threads with fingers. Add FGS-2 to the basal medium and mix well. Since several components are light-labile, the medium should not be exposed to light for extended periods. We do not recommend warming medium in a 37°C water bath prior to use. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.

#### Note

FGS-2 is for research use only. It is not approved for human or animal use, or for application in in vitro diagnostic procedures.